

 $f(x+\Delta x) = \sum_{i=0}^{\infty} \frac{(\Delta x)^i}{i!} f^{ij}(x) \int_a^b \int_a^b \frac{17}{12} \left\{ 2.71828 \right\}$ 

Center for Fast Ultrasound Imaging

1

### **CT** reconstruction - outline


- CT scanners
- Projection and Radon transform
- Projection demo
- Fourier slice theorem
- Inverse Radon transform filtered backprojection
- Selection of filters
- Filtered backprojection algorithm
- Reading material: Prince & Links chapter 6





Center for Fast Ultrasound Imaging, Department of Health Technology Technical University of Denmark


2


DTU











.

DTU

### **Radon transform**


$$p(x',\phi) = \int_{-\infty}^{+\infty} f(x'\cos\phi - y'\sin\phi, x'\sin\phi + y'\cos\phi)dy'$$

f(x,y) – Attenuation image


x', y' – Gantry coordinate system


x, y – Patient coordinate system


Center for Fast Ultrasound Imaging, Department of Health Technology Technical University of Denmark 8/x





C


























#### **Fourier Slice Theorem**

$$P(\rho,\phi) = \int_{-\infty}^{+\infty} p(x',\phi)e^{-j2\pi\rho x'} dx'$$
$$= F(\rho\cos\phi, \rho\sin\phi)$$

F(u,v) – 2D Fourier transform of image  $\phi$  – Gantry rotation

Center for Fast Ultrasound Imaging, Department of Health Technology Technical University of Denmark

20/x



## Filtered Back Projection (FBP)

$$\hat{f}(x,y) = \int_{0}^{\pi} \int_{-\infty}^{+\infty} |\rho| P(\rho,\phi) e^{j2\pi\rho x'} d\rho d\phi$$

 $\hat{f}(x, y)$  – Reconstructed image  $\phi$  – Gantry rotation x' – Detector position


Center for Fast Ultrasound Imaging, Department of Health Technology

21/x


21


## Filtered backprojection


- Perform for all projections:
  - 1. Make Fourier transform of projected data
  - 2. Apply filter in Fourier domain
  - 3. Make invers Fourier transform
  - 4. Backproject and sum with previous image





enter for Fast Ultrasound Imaging, Department of Health Technolog echnical University of Denmark
































# Reading DICOM data 1



```
% Set overall parameters
dir name='DICOM/ST00001/SE00001/'; % Directory name
      start image=1; % First image in series
      end_image=747;
                                           % Last image in series
                                           % Frame rate for playing back the movie
      frame rate=50;
     bone off set= -250; % Offset for showing the bones bone_range=100; % Range for showing the bones
      % Initialize figure
      colormap(gray(map_values));
      dicom_movie(end_image+1-start_image) = struct('cdata',[],'colormap',[]);
Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark
```

Center for Fast Ultrasound Imaging, Department of Health Technology Technical University of Denmark

# Reading DICOM data 2



```
Read information for the first images
file_name='IM00001';
info=dicominfo([dir name, file name]);
dx=info.PixelSpacing(1);
dy=info.PixelSpacing(2);
Y = dicomread(info);
[Nx, Ny] = size(Y);
% Make space for all the images
Y=zeros(Nx,Ny,end_image+1-start_image);
z_positions=zeros(end_image+1-start_image,1);
```

36

## Reading DICOM data 3



```
% Loop through the images and read and display them
      for i=start image:end image
          file name=['IM00',num2str(floor(i/100)),num2str(floor(rem(i,100)/10)),num2str(rem(i,10))];
           info=dicominfo([dir_name, file_name]);
          Y(:,:,i) = dicomread(info);
          image(((1:Nx)-Nx/2)*dx,((1:Ny)-Ny/2)*dy,(double(Y(:,:,i))+off_set)/range*map_values)
          xlabel('Distance [mm]')
          ylabel('Distance [mm]')
          pos=sprintf('%5.1f, %5.1f, %5.1f', info.ImagePositionPatient(1), ...
info.ImagePositionPatient(2), info.ImagePositionPatient(3));
          z positions(i)= info.ImagePositionPatient(3);
text(-150, -150, ['Image', num2str(i),', Pos. (x,y,z) = ', pos,' mm'], 'Color', [1 1 1])
          axis('image')
          drawnow
          dicom movie(i)=getframe;
      % Display the movie
     movie(dicom_movie, 5, frame_rate);
Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark
                                                                                                                   Script in: ct_data/dicom_data
```

37



Full script can be found at:

https://courses.healthtech.dtu.dk/22485/files/ct\_data/dicom\_data/display\_dicom\_images.m

on the page for the CT data: https://courses.healthtech.dtu.dk/22485/?ct\_data/assign\_data.html

Center for Fast Ultrasound Imaging, Department of Health Technology Technical University of Denmark

38/x